



JDENTAL CARE

just smile


JD BONE TRACK™ DRILLS KIT

For a perfect post extractive implant insertion



Place implants in post extractive sites correctly and accurately with JD Bone Track drills!

The JD Bone Track drill is a registered product, specifically designed to simplify the insertion of the post-extractive implant. JD Bone Track drills are characterized by a diamond cutting body and a non-cutting tip. Follow the Bone Track method!



The JD Bone Track drills are specifically designed to simplify the immediate insertion of the post-extractive implant.

JD Bone Track drills are characterized by a diamond cutting body and a non-cutting tip. They are available in four different implant diameters (Ø 3.2, Ø 3.7, Ø 4.3, Ø 5.0mm).

→ Diamond cutting body
→ Non-cutting tip

Note: all dimension are expressed in millimeters

Follow the Bone Track Method


A new approach that will simplify your daily clinical practice. Insert the non-cutting tip of the JD BoneTrack Drill into the underprepared osteotomy and push the drill palatally in order to create a track on the palatal bone. In this way, the correct space for the implant body will be created.

IMPLANT DIAMETER	PREPARATION IN MAXILLA	PREPARATION IN MANDIBLE
Ø 3,2	1.5 2.0 2.4 JD Bone Track Drill 3.2	1.5 2.0 2.4 2.8 JD Bone Track Drill 3.2
Ø 3,7	1.5 2.0 2.4 JD Bone Track Drill 3.7	1.5 2.0 2.4 2.4 2.8 JD Bone Track Drill 3.7
Ø 4,3	1.5 2.0 2.4 2.8 JD Bone Track Drill 4.3	1.5 2.0 2.4 2.4 2.8 3.2 JD Bone Track Drill 4.3
Ø 5,0	1.5 2.0 2.4 2.8 3.2 JD Bone Track Drill 5.0	1.5 2.0 2.0 2.4 2.8 2.8 3.2 3.6 JD Bone Track Drill 5.0





1. Identify the anatomy of the socket and start drilling towards the palatal bone.

N.B.: The last standard drill diameter to be used for under-prepare the osteotomy always corresponds to the implant tip diameter to be inserted.



2. Insert the non-cutting tip of the JD BoneTrack Drill into the osteotomy and push the drill palatally in order to create a track on the palatal bone.



3. Insert the implant leaving the correct vestibular gap

Bone Track Drills:

JDDIADR32	Diamond Drill Ø 3.2
JDDIADR37	Diamond Drill Ø 3.7
JDDIADR43	Diamond Drill Ø 4.3
JDDIADR50	Diamond Drill Ø 5.0



JD Bone Track Drills Kit:

JDBTK	JD Bone Track Drill Kit
-------	-------------------------



JDentalCare, the JDentalCare logo and other trademarks mentioned in this document are, excluding specifically mentioned exceptions, trademarks property of JDentalCare.

The images of the products represented in this brochure are not necessarily to scale.

Note: Not all products manufactured by JDentalCare and represented in this brochure have been regulatory cleared/released for sale in all countries and/or have been licensed in accordance with FDA Rules and Regulations.

For more information about current product assortment and availability, please contact your local JDentalCare distributor. For prescription only.

Caution: Federal (US) law restricts this device to sale by or on the other for a licensed dentist. See instructions for use for more information such as indications, warnings, precaution and contraindications.

REV.02 of 2024-05-08