

JD GUIDED SURGERY WORKFLOW

JDentalCare accompanies you all along the way and puts at your disposal the experience and capabilities of JDLab, a laboratory of the latest generation to which you can turn

1

THE DENTIST SENDS TO THE JDLAB:

- CT scan in DICOM format
- Diagnostic wax-up copy
- Templates
- Medical prescription
- Evobite

2

THE JDLAB PREPARES:

- Model scans with Evobite
- Guided-surgery digital planning
- Creation of the surgical guide

3

THE JDLAB SENDS TO DENTIST:

- CT scan in DICOM format
- Diagnostic wax-up copy
- Templates
- Medical Prescription
- Evobite
- +
- **Surgical guide**
- **Chewing index**
- **3D printed model**
- **Digital planning summary**

